



ET Tape™

**Minimizes inadvertent extubations;
ensures tube stability**

ET Tape is a pre-cut, non-irritating foam adhesive endotracheal tube holder designed to ensure stability of the endotracheal tube and minimize inadvertent extubations. ET Tape's latex free hypoallergenic materials have superior holding power and protect the upper lip and skin across the cheek area. A skin prep wipe is provided to form a topical barrier between skin and adhesive for added protection against skin irritation and breakdown. A Releasable Cable Tie secures placement of the endotracheal tube and prevents slippage. ET Tape provides tube stability for short-term applications and allows for consistency of taping and securement. ET Tape comes as a prepackaged kit that includes tape, skin prep and Releasable Cable Tie.

B&B recommends using Universal Bite Block™ or Bite Proof Bite Block™ with ET Tape to protect teeth and prevent biting on the endotracheal tube. ET Tape can be used in all critical care units, operating rooms, emergency departments, special procedure units and transport situations.

Ordering Information

11020	ET Tape (100/case)
11025	ET Tape (50/case)

Complementary Products

11160	Universal Bite Block (Adult) (50/case)
11170	Universal Bite Block (Pediatric) (50/case)
11040	Bite Proof Bite Block (100/case)

B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.
Visit www.BandB-Medical.com or contact us today at 1.800.242.8778.

SAFE

Latex free, hypoallergenic

Skin prep provides topical barrier to minimize skin irritation and breakdown

Enhanced durability for endotracheal tube protection and position

Non absorbent material reduces tube slippage

AHA guidelines recommend use of a commercial pre-cut device for ET tube security

COST EFFECTIVE

Prepackaged kit for immediate use includes pre-cut tape, skin prep and Releasable Cable Tie

Reduces cost of repeated taping procedures

CONVENIENT

Pre-cut tape minimizes application time

Easy and quick application for time sensitive conditions

VERSATILE

Used for endotracheal tubes 4 mm and greater

Child through adult application

Oral or nasal application

Designed for use in critical care and special procedure units, operating rooms, emergency and transport situations

Accommodates patients weighing 40 to 300 pounds

11020C2

© 2009 B&B Medical Technologies. All rights reserved.
Skin-Prep™ is a trademark of Smith & Nephew.

"A World of Products for Better Breathing" is a service mark of B&B Medical Technologies.



Toll-free: +1.800.242.8778

Tel: +1.760.929.9972

Fax: +1.760.929.9953

2734 Loker Avenue West, Suite M
Carlsbad, CA 92010 USA

