



# Universal Bite Block™

Prevents ET tube occlusion while biting;  
protects cuff inflating line

The Universal Bite Block is the most complete system to prevent biting of the endotracheal tube and occluding or rupturing the cuff inflating tube. The latex free Universal Bite Block secures directly on the endotracheal tube and allows for delivery of proper oral hygiene by providing unobstructed bilateral access to the oral cavity. The Universal Bite Block can also be used for protection of closed suction catheters, endoscopes, bronchoscopes and oral gastric tubes.



Universal Bite Block (Adult)

The adult-sized Universal Bite Block fits 6 mm and larger endotracheal tubes and is secured with a Releasable Cable Tie. The pediatric-sized Universal Bite Block fits 4 to 6 mm endotracheal tubes and is secured with a cable tie. When

utilized with B&B's StabiTube™, LockTite™ and ET Tapes™, or with other products and adhesives used to secure and support the artificial airway, the Universal Bite Block provides the fastest and most secure system for protecting the artificial airway.



Universal Bite Block (Pediatric)

## SAFE

Latex free

Hypoallergenic materials for fragile skin

Provides convenient access for oral care, suctioning and procedures

Well tolerated by patients

## COST EFFECTIVE

Choice for short and long term applications

Packaged with securing device

## CONVENIENT

Easy and quick application for time sensitive conditions

Use with B&B and all endotracheal tube holders

## VERSATILE

Multiple sizes fit 4 to 6 mm tubes and 6 mm and larger tubes

Designed for use in critical care, surgery, emergency, transport and alternate care environments

Can be used in combination with B&B endotracheal tube holders and other commercially available securing devices

## Ordering Information

11160 Universal Bite Block Adult (50/case)

11170 Universal Bite Block Pediatric (50/case)

## Complementary Products

11060 StabiTube (40/case)

11010 ET Tape II (100/case)

11020 ET Tape (100/case)

11090 LockTite (50/case)

## B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.

Visit [www.BandB-Medical.com](http://www.BandB-Medical.com) or contact us today at 1.800.242.8778.

11160C2

© 2009 B&B Medical Technologies. All rights reserved.

"A World of Products for Better Breathing" is a service mark of B&B Medical Technologies.

Toll-free: +1.800.242.8778

Tel: +1.760.929.9972

Fax: +1.760.929.9953

2734 Loker Avenue West, Suite M  
Carlsbad, CA 92010 USA

