



TrachGuard™

Complete kit prevents inadvertent disconnects; stabilizes ventilator circuit and tracheostomy tube

TrachGuard is a patented combination long term tracheostomy tube holder and anti-disconnect device that dramatically reduces the potential for inadvertent ventilator circuit disconnects and exposure to cross contamination. Made to accommodate any tracheostomy tube, the anti-disconnect device stabilizes the patient connection to the ventilator circuit or in-line suction device.

TrachGuard is a complete kit that secures the tracheostomy tube to the patient connector of the ventilator circuit and a comfortable soft collar with hook and loop closure to secure the tracheostomy tube in position. The TrachGuard collar has superior holding power and is made of latex free, hypoallergenic soft foam. The collar fits neck sizes 7 to 18 inches. TrachGuard is an essential tool for busy ICUs, long term acute care, home and alternate care facilities.



SAFE

Latex free, hypoallergenic materials

Enhanced durability for securing the ventilator circuit connection

Minimizes potential for decannulation

Reduces potential exposure to contaminated tubing moisture

COST EFFECTIVE

Prepackaged kit with collar and anti-disconnect device

Designed for use in critical care, PICU, long term acute care, alternate site and home care

CONVENIENT

Prepackaged kit minimizes time for setup and connection

Easy and quick application to all tracheostomy tubes

Comfortable foam collar fits pediatric through adult patients

VERSATILE

Adult and pediatric applications

Collar fits neck sizes 7 to 18 inches

Can be used with delivery of aerosol and supplemental oxygen therapy

Ordering Information

11250 TrachGuard (80/case)

11255 TrachGuard (20/box)

Complementary Products

20001 Sil.Flex Stoma Pad™ – Size 1 (10/box)

20002 Sil.Flex Stoma Pad – Size 2 (10/box)

20003 Sil.Flex Stoma Pad – Size 3 (10/box)

B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.

Visit www.BandB-Medical.com or contact us today at 1.800.242.8778.



11250C2

© 2009 B&B Medical Technologies. All rights reserved.

Sil.Flex Stoma Pad is a trademark of A Plus Medical.

"A World of Products for Better Breathing" is a service mark of B&B Medical Technologies.

Toll-free: +1.800.242.8778

Tel: +1.760.929.9972

Fax: +1.760.929.9953

2734 Loker Avenue West, Suite M
Carlsbad, CA 92010 USA

