



# HOPE™ Nebulizer

**Delivers outstanding nebulizer performance with patented Heliox port**

The HOPE Nebulizer has proven to be a cost effective solution for delivery of continuous medication for the treatment of moderate to severe asthma and compromised airways. The HOPE is a high output nebulizer used for delivery of aerosolized medications and diagnostic formulations. It provides additional hydration for patients to loosen secretions when extended therapy is required.

The HOPE is a refillable, closed dilution, high output, large volume continuous medication nebulizer. The flow design and 200 ml reservoir allow for up to 8 hours of nebulized medication delivery at 25 ml/hr. The HOPE delivers 1 to 8 hours of continuous medication with 2.5 MMAD. The output at 25 ml/hr delivers >60% particles in the respirable range at 10 lpm flow. With the addition of supplemental Heliox, the HOPE provides a 50 ml/hr output with >75% particles in the respirable range with 15 lpm supplemental 80/20 Heliox delivered through the auxiliary port.

The HOPE is the gold standard for delivery of continuous nebulization and Heliox therapy. Heliox has been reported to build a therapeutic bridge in severe asthma and compromised airways until conventional medications have had time to take effect. B&B's patented supplemental port for Heliox delivery greatly enhances the performance of the HOPE Nebulizer. It is no wonder that the competition has reconfigured their products to emulate the HOPE Nebulizer. B&B first introduced the patented supplemental port for Heliox delivery in a closed dilution nebulizer, greatly enhancing the performance of the HOPE Nebulizer. HOPE is the industry leader in simplicity and outstanding nebulizer performance in a Large Volume Continuous Medication Nebulizer.

## SAFE

Rescue therapy for asthmatics and obstructive lung disease patients

Latex-free solution

Easy-to-understand dosing chart

Allows clinician to add auxiliary gas flow without recalculation of medication

Hypoallergenic materials

## COST EFFECTIVE

Prepackaged kit with needed accessories for delivery of continuous medication therapy

Dosing chart provided with each HOPE Nebulizer

Delivers medication therapy and humidification for up to 8 hours

## CONVENIENT

Quick connection to flowmeter

Easy and quick application for time-sensitive conditions

Twist off reservoir for easy medication changes

Standard auxiliary port for delivery of therapeutic gases

## VERSATILE

Designed for use in the critical care units, emergency department, extended care clinics, outpatient clinics and in hospital transport

Accommodates auxiliary port for delivery of supplemental heliox gas for therapy

## Ordering Information

11310	<b>HOPE Nebulizer</b> (25/case)
11320	<b>HOPE Nebulizer Adult Kit</b> (10/case) Nebulizer, 6' Aerosol Tubing and Adult Mask
11330	<b>HOPE Nebulizer Pediatric Kit</b> (10/case) Nebulizer, 6' Aerosol Tubing and Pediatric Mask
11370	<b>Heliox Regulator</b> (1 each) With Flowmeter and Nipple
11373	<b>70/30 Heliox DISS Regulator</b> (1 each) With Flowmeter, Nipple and DISS check valve
11382	<b>80/20 Heliox DISS Regulator</b> (1 each) With Flowmeter, Nipple and DISS check valve



## B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.

Visit [www.BandB-Medical.com](http://www.BandB-Medical.com) or contact us today at 1.800.242.8778.

Toll-free: +1.800.242.8778

Tel: +1.760.929.9972

Fax: +1.760.929.9953

2734 Loker Avenue West, Suite M  
Carlsbad, CA 92010 USA



11310C4

© 2011 B&B Medical Technologies. All rights reserved. The HOPE Nebulizer is a trademark and "A World of Products for Better Breathing" is a service mark of B&B Medical Technologies.