



The Test Lung 1.0 L (Adult)

## The Test Lung™

**Simulates nominal compliance and resistance;  
included connectors adapt to all patient circuits**

The Test Lung 1.0 L (Adult) and Test Lung 0.5 L (Pediatric) provide an economical choice for providing high quality demonstration and testing applications on mechanical ventilators for adult or pediatric patients. The Test Lung simulates the respiratory system, providing nominal levels of resistance and compliance. The connector kit includes a variable leak adjustment to demonstrate patient-trigger function or to simulate an airway leak. The Test Lungs are packaged with a 1.0 or 0.5 liter, latex free silicone ventilation bag, Test Lung Connector Kit™ and Rp5 Resistance. Each bag is durable, easily removable and can be cleaned or sterilized as needed. The Test Lung Connector Kit comes standard with each Test Lung and adapts to all patient circuits and proximal airway flow sensors. The kit contains one elbow connector and two straight connectors with leur ports and locks.

The Test Lung is compact in design and lightweight. Each Test Lung is validated for resistance and compliance in the application range and has a unique serial number to ensure its compliance with specification. The Test Lung is the ideal tool for teaching and demonstration as well as for performing ventilator verification testing in respiratory care, biomedical labs and anesthesia departments.

### Ordering Information

20110 The Test Lung 1.0 L (Adult) (1/box)

20120 The Test Lung 0.5 L (Pediatric) (1/box)

### Complementary Product

20118 Precision Resistor Kit™ (1/box)

### B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.

Visit [www.BandB-Medical.com](http://www.BandB-Medical.com) or contact us today at 1.800.242.8778.



The Test Lung 0.5 L (Pediatric)



### SAFE

Latex free, hypoallergenic

Factory calibrated

Can be cleaned and sterilized between uses

### COST EFFECTIVE

Complete kit for demonstration of ventilator, anesthesia equipment and circuit leaks

Silicone ventilation bag easily removed for cleaning

### CONVENIENT

Packaged with adapters for ventilators

Easy to use, no assembly required

### VERSATILE

1.0 L silicone bag for adult;  
0.5 L silicone bag for pediatric

Designed for use with demonstration and testing of ventilators in respiratory care, clinical teaching areas and biomedical departments

20110C2

© 2009 B&B Medical Technologies. All rights reserved.

"A World of Products for Better Breathing" is a service mark of B&B Medical Technologies.

Toll-free: +1.800.242.8778

Tel: +1.760.929.9972

Fax: +1.760.929.9953

2734 Loker Avenue West, Suite M  
Carlsbad, CA 92010 USA

