



Danny Ties™

Unique tracheostomy holder offers comfortable fit; quilted design minimizes skin irritation under the collar

To ensure safety, the patient with a tracheostomy needs to have a tube holder that securely holds the artificial airway in place so that the device is not accidentally dislodged. As important, the holder needs to provide a soft and comfortable fit about the neck without irritating the skin under the collar. Danny Ties are made of a soft, absorbent material that lays smoothly at the edges of the collar and significantly minimizes skin breakdown beneath the collar. The collar holds its shape, does not fold in half around the neck and does not stretch when it absorbs moisture. The patent-pending design of the Danny Ties evenly distributes the padded collar material around the neck to minimize pressure points on the skin.

The Danny Ties collar is easy to apply with tapered ends to the collar straps. The ends thread easily through the eyelets of the tracheostomy tube, allowing for quick application and changes of the collar on the smallest infants and the largest adult patients. If soiled, Danny Ties may be washed, dried and reused for the same patient per institutional policy.

Danny Ties were developed with care by a dad wanting to “make a difference” for his son and change the quality of life for tracheostomy patients. His engineering expertise and commitment to finding a better solution is found in the Danny Ties.

The Sil.Flex™ Stoma Pad may be used in combination with Danny Ties. Sil.Flex Stoma Pads enhance patient comfort with all brands and styles of tracheostomy tubes. The TrachStay™ may also be used in tandem with the Danny Ties to prevent accidental disconnects from the ventilator circuit.

Ordering Information

Small (fits neck 6 to 10 inches)

25001 Danny Ties – Small (5/bag)

25011 Danny Ties – Small (10/box)

25101 Danny Ties – Small (100/cs)

Large (fits neck 19 to 25 inches)

25003 Danny Ties – Large (5/bag)

25013 Danny Ties – Large (10/box)

25103 Danny Ties – Large (100/cs)

Medium (fits neck 9 to 20 inches)

25002 Danny Ties – Medium (5/bag)

25012 Danny Ties – Medium (10/box)

25102 Danny Ties – Medium (100/cs)



B&B Products are available from

B&B Medical Technologies and finer specialty care distributors.

Visit www.BandB-Medical.com or contact us today at 1.800.242.8778.

SAFE

Quilted collar provides a softer, more comfortable fit about the neck

Absorbent collar helps reduce skin breakdown around the neck

Environmentally friendly packaging

Latex free, hypoallergenic tracheostomy collars

COST EFFECTIVE

28 day application, Danny Ties may be reprocessed per institutional policy

Packaged 5 Danny Ties per pack

CONVENIENT

Unique strap design allows for easy application and removal of the tie

Single Patient Use

Available in three sizes to accommodate infant through adult application

VERSATILE

Tapered straps for easy application through tracheostomy tube eyelets

Unique sizes and contoured shape of the fastener to fit different neck sizes

Special design minimizes rolling and stretching of the collar



25000C3

© 2012 B&B Medical Technologies. All rights reserved.
 Danny Ties is a trademark of Leckie Medical Products, Inc.
 “A World of Products for Better Breathing” is a service mark of B&B Medical Technologies.

Toll-free: +1.800.242.8778

Tel: +1.760.929.9972

Fax: +1.760.929.9953

1954 Kellogg Avenue
 Carlsbad, CA 92008 USA

